



FCI®

Your Patients, Our Expertise

PRECIVIA®

**THE GOLD STANDARD*
FOR ALL INTRAVITREAL INJECTIONS**



**FCI
EXCLUSIVE**

*Petri, A., Boysen, K., Cehofski, L.J. et al. Intravitreal Injections with Vascular Endothelial Growth Factor Inhibitors: A Practical Approach. Ophthalmol Ther 9, 191-203 (2020).

**MORE THAN 2 MILLION INJECTIONS
PERFORMED WORLDWIDE!**



PRECIVIA®

Single use

Packaging

S5.200.P

To be used with syringes filled with solutions and needles 30 Ga x 1/2 inches, connected by a luer-lock connection.

Box of 25

S5.200.I

Compatible with insulin syringes with crimped needle.

Box of 25

PRECIVIA® PACK

- 1x Wrapping Sheet 75x75 cm
- 1x Perforated Drape 53x50 cm
- 1x S5.200.P Precivia®
- 3x Pro Ophtha Sticks
- 5x Compresses 5x5 cm
- 1x 30Ga 12 mm Cannula



Single use

Packaging

S5.300.P

Pack IVT 1

Box of 25

Precivia® is a Class Is medical device indicated for intravitreal injections and manufactured by FCI S.A.S. - Notified Body: GMED CE n°0459.
Precivia® Pack is a Class IIa medical device indicated for intravitreal injections. It is distributed by FCI S.A.S. and manufactured by M.E.D. Medical Products GmbH - Notified Body: DQ5 Medizinprodukte GmbH CE n°0297.

FCI20/006 ind.B - October 2021 - Société par Actions Simplifiées au Capital de 49909 € - SIREN 351 670 385 - RCS Paris.



fciworldwide.com



fciworldwide.com/contact